04a AP Determination of % Fe by Redox Titration (7951989)

Question

1 2 3 4 5 6 7 8 9

Description

Basic Periodic Table

Instructions

Determination of % Fe by Redox Titration

Introduction: Acid-base titration is used to determine the quantitative relationship between two reacting solutions containing hydrogen ions and hydroxide ions. Titration methods can also be used to determine the amount of other ions in solution if the substance is oxidized or reduced to give a color change. Such a method works well to determine the amount of iron in an iron(II) salt using permanganate because the permanganate ion (MnO_4^-) is purple, but the Mn^{2+} is colorless: (unbalanced)

 $Fe^{2+} + MnO_4^-$ (purple) --> $Fe^{3+} + Mn^{2+}$ (colorless)

Materials: 250 mL flask, repipet, buret, balance, magnifying glass, index card, unknown Fe^{2+} salt, 0.0200 M KMnO₄, 6.0 N H₂SO₄

Objective: To determine the % Fe in an unknown salt.

Procedure:

1. Weigh out between 0.500 - 0.600 g of the unknown iron salt on the electronic balance using a plastic weighing dish. Record the mass. Transfer the salt to a 250 mL erlenmeyer flask.

2. Dissolve the iron salt in about 50 mL of distilled water in the 250 mL erlenmeyer flask. (use the approximate graduations on the flask to measure the 50 mL of water) 3. While wearing goggles, using the repipet, obtain10.0 mL of 6.0 N H_2SO_4 in the flask

with the iron salt. Swirl to mix.

4. Titrate with 0.0200 M KMnO₄ until a faint pink color persists for 30 seconds. (if a brown precipitate forms, start over and make sure the acid was added) Discard the solution in the sink.

5. Repeat steps 1-4 for 1 careful trial.

6. Rinse the buret with water and leave it upside down with the valve open.

7. Record your data and calculations in the computer.

```
Data:
```

Mass of unknown salt: Rapid Trial_____ Trial 1_____ Volume of KMnO₄: Final Reading Final Reading

Initial Reading______ Initial Reading______

Rapid Trial _____ Trial 1_____

Questions: (The questions below are in this sequence for a reason.)

1. Write the balanced equation for this reaction. The skeleton equation is written in the introduction. (Acid or base solution?)

2. Calculate the moles of MnO₄- using trial 1. (MnO₄- is 0.0200 M, K⁺ is a spectator ion)

3. Calculate the moles of iron in the salt sample. (Hint: you need the balanced equation)

- 4. Calculate the mass of Fe in the sample.
- 5. Calculate the % of Fe in the sample using trial 1.

1.	Question Details	Lab Partners [1837468]
	Enter the name(s) of your lab partner(s). (If	you worked by yourself, enter "none").
2.	Question Details	Objective and procedure summary [3413760]
	Restate the objective in your own words usin include any safety concerns).	g complete sentences. Summarize the steps in your procedure. (Be sure and
3.	Question Details Upload Lab Photo [3413757] Upload a photo of the lab apparatus with your face in the photo as you perform some part of the lab. Title the image with a unique file name before you upload it.(Maybe use your initials and part of the lab title) Choose File no file selected It must be less than 5 MB in size.	
4.	Question Details Write the balanced equation for the reaction your answer.)	Determination of %Fe by Redox Titration Balance [1743296]

5.	Question Details	Determination of %Fe by Redox Titration [1743290]	
	a. Enter the experimental mass(from the balance) of the unknown salt for trial 1: 40		
	b. Enter the initial reading of the buret for trial 1:40		
	c. Enter the final reading of the buret for trail $1:40$ mL		
	d. Calculate the volume of KMnO ₄ used in trial 1:		
	e. Calculate the moles of MnO ₄ ⁻ reacted in trial 1:49 mol		
	g. Calculate the mass of Fe in the sample:		
	h. Calculate the % of Fe in the sample:		
6.	Question Details	Percent Error % Fe by Redox Lab [1743603] _	
	Enter the percent of Fe in the unknown salt again (same number as question 5h):		
	You will probably get the red X on this question. The accuracy of your a	answer determines the number of points deducted.	
7.	Question Details	Observations, Skills utilized and learning [3413764]	
	What observations did you make during the lab? What chemistry concepts, laws, and/or skills were necessary to complete this		
	lab? What did you learn or re-learn? Use complete sentences.		
8.	Question Details	Error discussion [3413763]	
	What are some specific sources of error, and how do they influence the data? Which measurement was the least precise? Does the error make the final value obtained larger or smaller than it should be (give at least one example and trace the steps)? If your calculated percent errors are significant, you must propose valid explanations here. Instrumental error and human error exist in all experiments, and should not be mentioned as a source of error unless they caused a significant fault. Significant digits and mistakes in calculations are NOT a valid source of error. In writing this section		
	it is sometimes helpful to ask yourself what you would do differently if you were to repeat the experiment and wanted to obtain better precision and accuracy. Use complete sentences		

9. Question Details

Upload Calculations (Show Work) [3418656]

Upload a photo of your calculations, showing your work. Make sure your name and the date are written on the page. Title the image with a unique file name before you upload it.(Maybe use your initials and part of the lab title and the word Calcs) Choose File no file selected

Assignment Details

Name (AID): 04a AP Determination of % Fe by Redox Titration (7951989) Feedback Settings Submissions Allowed: 5 Before due date Category: Homework **Question Score** Code: Locked: Yes Author: Ryan, Matt (mryan@allsaintsschool.org) Last Saved: Sep 25, 2016 09:20 PM CDT Mark Permission: Protected Randomization: Person Help/Hints Which graded: Last Response Save Work After due date

Feedback Settings Before due date Question Score Assignment Score Publish Essay Scores Question Part Score Mark Add Practice Button Help/Hints Response Save Work After due date Question Score Assignment Score Publish Essay Scores Key Question Part Score Solution Mark Add Practice Button Help/Hints Response